


UVP ColonyDoc-It

Accurate and automated colony counting is possible with this compact station. High resolution megapixel digital color camera captures white light and fluorescent marked colonies. The epi blue light and the optional GFP filter enable users to visualize GFP fluorescence detection.



Dimensions: 13W x 12.5D x 17.5H in. (343 x 318 x 445 mm), 20 lbs.



For current camera and system specifications visit our website.

Order Number		Description
115 V	230 V	UVP ColonyDoc-It Imaging Station and Accessories
97-0539-01	97-0539-02	ColonyDoc-It Imaging Station
38-0340-01		Filter, GFP

Thermal Printer

Document gel and blot images with a 256 grayscale thermal printer. The compact printer easily connects to any of Analytik Jena's imaging systems (except the PhotoDoc-It).

- Low cost, prints in seconds
- High resolution, 3 x 4" archive quality thermal prints



Dimensions: 6W x 9.5D x 3.3H in. (154 x 240 x 84 mm), 7.5 lbs.

Order Number		Description
115 V	230 V	Thermal Printers and Accessories
89-0069-06	89-0069-07	Digital Thermal Printer, Mitsubishi
89-0038-01		Thermal Paper, Mitsubishi, High density, 4 rolls (800 prints)
89-0174-01		Thermal Paper, Sony, Super Glossy, UPP110HG 5 rolls (1000 prints)