

Portable & Rechargeable UV Lamps



Portable UV Lamps are lightweight for field use.

- **Portable UV Lamps** (6-watt) are operated by two 6-volt batteries
- The distinctive wrap around handle comfortably balances the lamp in your hand
- A built-in flashlight is included
- Optional rechargeable 6V batteries and a charger are available

Lamp Dimensions:
Portable: 9.5L x 2.8W x 9.4H in. (241 x 71 x 239 mm), 4 lbs.
Rechargeable: 5.5L x 3.0W x 9.8H in. (140 x 76 x 249 mm), 4 lbs.



Rechargeable UV Lamps with built-in rechargeable battery. UVSL-26P (shown) can operate with both longwave and shortwave at the same time.

- **Rechargeable UV Lamps** (6-watt) feature a lightweight design and come in three dual-tube models
- Lamps are equipped with an internal gel based, lead-acid battery which can be charged from a standard wall outlet
- Optional 12V adapter powers the lamp from a car cigarette lighter
- Use with the C-10P Cabinet (see page 35 for details)

| Order Number | Description |
|--|--|
| Portable UV Lamps, Replacement Tubes and Accessories | |
| 95-0011-01 (12V) | Portable UV Lamp, ML-49 6-watt, 365 nm BLB UV |
| 95-0009-02 (12V) | Portable UV Lamp, UVG-47 6-watt, 254 nm UV |
| 95-0010-02 (12V) | Portable UV Lamp, UVGL-48 split-tube, 6-watt, 254 nm/365 nm UV |
| 34-0016-01 | Replacement Tube, ML-49 |
| 34-0013-01 | Replacement Tube, UVG-47 |
| 34-0015-01 | Replacement Tube, UVGL-48 |
| 34-0014-01 | Bulb, White |
| 45-0005-01 | Battery, 6V |
| 45-0005-02 | Battery, rechargeable 6V (115V) |
| 58-0136-01 | Charger (115V) |

| Order Number | | Description |
|--------------|------------|--|
| 115 V | 230 V | Rechargeable UV Lamps, Replacement Tubes and Accessories |
| 95-0186-01 | 95-0186-02 | Rechargeable UV Lamp, UVL-26P, 6-watt, 365 nm UV |
| 95-0187-01 | 95-0187-02 | Rechargeable UV Lamp, UVS-26P 6-watt, 254 nm UV |
| 95-0181-01 | 95-0181-02 | Rechargeable UV Lamp, UVSL-26P 6-watt, 254 nm/365 nm UV |
| 34-0016-01 | | Replacement Tube, 365 nm BLB, UVL-26P |
| 34-0013-01 | | Replacement Tube, 254 nm, UVS-26P and UVSL-26P |
| 34-0034-01 | | Replacement Tube, 365 nm, UVSL-26P |
| 58-0127-02 | | Rechargeable UV Lamp Adapter |

Signature UV Lamps

Signature and fraud detection UV lamps use longwave UV to check the authenticity of signatures, documents, tickets, passports, currency and credit cards.

- Easy to operate
- Easy to detect forged documents and currency
- Safe to use

SL-2M Lamp: Invisible inks or altered documents fluoresce when exposed to the UV.



Dimensions:
SL-2M: 6.5W x 5.75D x 4.5H in. (165 x 145 x 114 mm), 3 lbs.

Mineralight UV Lamps

Mineralight R-52G Grid Lamp provides a powerful, uniform source of high-intensity 254 nm shortwave UV.

- For small scale mineral displays and laboratory applications including bacteria destruction and photochemistry
- Removable filter assembly converts the lamp for each of your application purposes

The 2UV models allow users to switch between 302 nm and 365 nm UV wavelengths in one unit.



Lamp Dimensions:
7.8H x 16L x 7.8W in. (198 x 406 x 198 mm);
Primary cord: 8 ft., secondary cord 3 ft.; Weight: 9.5 lb (4.3 kg)

| Order Number | | Description |
|--------------------------------|------------|--|
| 115 V | 230 V | Signature UV Lamps, Mineralight UV Lamps and Replacement Tubes |
| 95-0164-01 | 95-0164-02 | SL-2M Fraud Detection, UV Lamp, 4-watt, 365 nm UV |
| 95-0125-01 (4-AA not included) | | UVL-4 UV Lamp, 4-watt, 365 nm UV |
| 94-0012-05 | 95-0012-06 | R-52G Grid Lamp, 254 nm UV |
| 34-0010-01 | | Replacement Tube, 4-watt, BLB UV (UVL-4 Lamp) |
| 77-0002-02 | | R-52G Grid Replacement |