


Ultraviolet Lamps

Analytik Jena manufactures a large selection of Mineralight, Blak-Ray and UV lamps to accommodate many applications and budgets.

Ultraviolet Lamps: Hand & Specialty

Compact UV Lamps

Compact UV Lamps (4-watts) feature models with longwave (365 nm), shortwave (254 nm), or multiple UV. J-124 Stand for use with Compact UV Lamps.



Dimensions: 7.8L x 2.8W x 2.1H in. (198 x 71 x 53 mm), 1 lb.
J-124 Stand: 6.5H x 7.5W x 4.8D in. (165 x 190 x 122 mm), 1 lb



- Combination wavelength version comes with a snap-on selector plate for switching between wavelengths
- Lightweight, easy-to-handle lamp design
- Comes in 12V style (for cigarette lighter operation)

Order Number

Description

115 V	230 V	Compact UV Lamps, Replacement Tubes and Accessories
95-0016-14	95-0016-15	UVG-11, 4-watt, 254 nm UV
95-0017-09	95-0017-10	UVGL-15 (split-tube), 4-watt, 365 nm/254 nm UV
95-0018-02	95-0018-03	UVL-21, 4-watt, 365 nm UV
95-0021-12	95-0021-10	UVGL-25, 4-watt, 365 nm/254 nm UV
95-0019-01	95-0019-03	UVL-23, 4-watt, 365 nm BLB UV
95-0019-14 (12V)		UVL-23R, 4-watt, 365 nm BL UV
34-0003-01		Tube, 4-watt, 254 nm UV
34-0004-01		Tube, 4-watt, 365/254 nm UV
34-0005-01		Tube, 4-watt, 365 nm UV
34-0010-10		Tube, 4-watt, 365 nm, BLB UV
98-0020-03		J-124 Lamp Stand