

Handheld UV Lamps

Handheld UV Lamps (6-watts) include shortwave (254 nm), longwave (365 nm) UV or the combination shortwave/ longwave in one lamp. Pair with J-129 Lamp Stand for hands-free lamp use.

- Comfortable, ergonomically designed handle for easy handling
- On/off switch is located on the back of the lamp



Dimensions: 14.9L x 3.2W x 2.5H in., 2 lbs.  
J-129 Stand: 9.75H x 13L in. (248 x 330 mm), 3 lbs.

Order Number		Description
115 V	230 V	Handheld UV Lamps, Replacement Tubes, and Accessories
95-0004-09	95-0004-10	UVG-54, 6-watt, 254 nm UV
95-0005-05	95-0005-06	UVGL-55 (split-tube), 6-watt, 365 nm/254 nm UV
95-0006-02	95-0006-03	UVL-56, 6-watt, 365 nm UV
95-0007-05	95-0007-06	UVGL-58, 6-watt, 365 nm/254 nm UV
95-0104-01	95-0104-02	UVM-57, 6-watt, 302 nm UV
34-0013-01		Tube, 6-watt, 254 nm UV
34-0015-01		Tube, 6-watt, 365/254 nm UV
34-0034-01		Tube, 6-watt, 365 nm UV
34-0044-01		Tube, 6-watt, 302 nm UV
98-0016-03		J-129 Lamp Stand

3UV Multi-Wavelength UV Lamps

This innovative 3UV lamp (patented) houses three UV wavelengths in one lamp. 3UV Lamps hold longwave, shortwave and midrange UV in one lamp (8-watt model pictured).



- Selection of each wavelength is by an easy turn of the dial located on one end of the housing
- The uniquely designed reflectors mounted behind each tube provide maximum UV for fluorescence applications
- Lightweight plastic and ergonomically-designed housing for handheld or stationary use

Lamp Dimensions:  
3UV-34: 9.5L x 3W x 4.5D in. (241 x76 x114 mm), 5 lbs.  
3UV-36: 12.5L x 3W x 4.5D in. (318 x76 x114 mm), 5 lbs.  
3UV-38: 15.5L x 3W x 4.5D in. (394 x76 x114 mm), 6 lbs.