

White/UV and Twin EL Series Lamps



UV/White and Twin EL Series Lamps are available in four, six and eight watt configurations. Lamps include UV and white light (UVS-14 shown).

- White/UV Models include one UV tube (254, 302 or 365 nm) and one white light tube in one lamp
- Twin Models include two tubes of the same UV wavelength for added UV intensity



Accessories for EL Series UV Lamps



J-138 Dimensions: 10.74H x 13L x 3W in. (273 x 330 x 76 mm)

EL Series Twin Models have two tubes of the same UV wavelength (UVL-28 shown).

- **J-138 Lamp Stand** is for use with any EL Series Lamps; allows for hands free use
- **C-10E4 and C-10E6 Cabinets** can be used with EL Lamps which emit UV in a darkroom environment (see page 35 for cabinet ordering information)
- The uniquely designed reflectors mounted behind each tube provide maximum UV for fluorescence applications
- The 8-watt EL lamps can be used with the **C-75 Cabinet** when a darkroom environment is needed (see pages 34-35 for cabinet ordering information)

Order Number		Description
115 V	230 V	UV/White EL Series Lamps and Lamp Accessories
95-0198-01	95-0198-02	UVL-18, 8-watt, 365 nm UV/white
95-0199-01	95-0199-02	UVM-18, 8-watt, 302 nm UV/white
95-0200-01	95-0200-02	UVS-18, 8-watt, 254 nm UV/white
95-0272-01	95-0272-02	UVL-16, 6-watt, 365 nm UV/white
95-0273-01	95-0273-02	UVM-16, 6-watt, 302 nm UV/white
95-0274-01	95-0274-02	UVS-16, 6-watt, 254 nm UV/white
95-0264-01	95-0264-02	UVL-14, 4-watt, 365 nm UV/white
95-0266-01	95-0266-02	UVS-14, 4-watt, 254 nm UV/white
18-0063-01		J-138 Lamp Stand
95-0072-08		C-10E4 Cabinet, (for 4W EL Lamps)
95-0072-09		C-10E6 Cabinet, (for 6W EL Lamps)